Experiment	Microbiota of the nose
Advisors	Prof. Dr. Brigitte Berger-Bächi, <u>bberger@immv.unizh.ch</u> , Department of Medical Microbiology, University of Zürich, Gloriastrasse 32, 8006 Zürich, Tel 044 634'26'50 Fax 044 634'49'06
Reading	Chapters in BBOM 10th: BBOM : Madigan M.T., J.M. Martinko and J. Parker: "Brock - Biology of Microorganisms", 10 <sup>th</sup> Edition (2003), Prentice Hall. <b>21.1-21.8</b> , 21.9-21.14, <b>26.9</b>
Objectives	To analyze the normal flora of the nose on selective media and identify staphylococci by biochemical tests.
Background	Normal bacterial flora. Body surfaces such as skin and mucous membranes are colonized by large numbers of microorganisms that are highly adapted to their particular environment. These organisms, also called commensals, constitute the body's normal flora. Positive effects of commensals are the continuous stimulation of the immune system of the host and protection from pathogens. When members of the normal flora from one anatomical site gain access to an otherwise sterile location, then infection may result. Sometimes these are relatively trivial, but others may be life-threatening. Such infections are known as endogenous. Predisposition for endogenous infections is destruction of the natural barriers against infection, such as lesions of the skin, or reduced immunity, or carriage of indwelling devices. Figure 1. Some bacteria of the skin and upper respiratory tract. Characteristic shape and staining properties. Blue, Gram-positive; red, Gram-negative. Skin
	<ul> <li>Coagulase-negative staphylococci</li> <li>Proionibacteria</li> </ul>
	• Corynebacteria
	• Yeast, fungi, especially <i>Candida albicans</i>





Department of Medical Microbiology, University of Zürich, Gloriastr.. 30/32, CH-8006 Zürich, URL: <u>http://www.imm.unizh.ch</u> Fax. 044 63 44906 Tel 044 63 42650

Practical Work	First week:						
I I ACULAI VVUIK	Each student wi	ill obtain:					
	a sheep	a sheep-blood plate					
	a CNA	plate					
	a McCo	onkey plate					
	Write your nam	Write your name on the back of each plate.					
	Take a sterile c into your nose cotton tip gently 3 <sup>rd</sup> fraction. Inc	Take a sterile cotton tip, dip it into sterile NaCl solution, and insert it about 1-2 cm into your nose and rub the interior nare with a light turning movement. Roll the cotton tip gently over about 1/3 of the plate (1 <sup>st</sup> fraction). Fractionate it into a 2 <sup>nd</sup> and 3 <sup>rd</sup> fraction. Incubate the plates at 37 °C. <b>Next week:</b> describe the colony morphology color amount of bacteria					
	Next week:						
	describe the col						
	color						
	amount of bacte						
	hemolysis type						
		Morphology of the colonies	Hemolysis				
	sheep blood						
	CNA						
	McConkey						
	Biochemical characterisation of the isolates.         Fach group will receive one S. aureus and one Staphylococcus epidermidis strain.						
	you have own	you have own colonies from the nose that look like staphylococci and show $\alpha$ -					
	hemolysis on th	hemolysis on the sheep blood plate, analyse them as well.					
	Diagram for th	Diagram for the identification of Gram-positive cocci					
		Gram-positive cocci					
		in clusters catalase +	in chains catalase -				
		Staphylococcus	Streptococcus				
	1						
	coagula	ase + coagulase -	γ-				
	$\alpha$ -hemo	biysis no hemolysis $\alpha$ -hen	noiysis β-nemolysis no hemolysis ↓				
	S. aure	us coagulase-neg. staphylococci	S. pyogenes (group A antigen) S. agalactiae (group B antigen)				

Department of Medical Microbiology, University of Zürich, Gloriastr.. 30/32, CH-8006 Zürich, URL: <u>http://www.imm.unizh.ch</u> Fax. 044 63 44906 Tel 044 63 42650

Catalase The catalase degrade	s the toxic hydrogen	peroxide into water a	and oxygen.	
	$2 H_2O_2 \rightarrow$	$\sim 2 \text{ H}_2\text{O} + \text{O}_2$		
The catalase is a defe by certain immune ce	ence mechanism of ba ells of the host to kill	acteria against the pe bacterial invaders (s	eroxide that is produce ee BBOM Fig. 6.29)	ed
Mix a loop full of ba appearance of bubble	ecteria with a drop of es is indicative of the	f 3 % H <sub>2</sub> O <sub>2</sub> on a glas presence of catalase	s slide. The immedia	ıte
<b>Tube coagulase test</b> The coagulase is an e	enzyme that causes cl	lotting of blood plasr	na.	
Mix a loop full of ba for fibrin clotting. Th the medium, while co	acteria with plasma, i ne pathogenic <i>S. auro</i> pagulase-negative sta	ncubate for several h eus produces coagula phylococci produce	ours at 37 °C and loo ase that is excreted in no coagulase.	ok to
<b>3c. Agglutination te</b>	st for <i>S. aureus</i> (Sta	phaurex Plus):		
The surface-bound agglutination of the agglutination kit reco factor and staphyloco factor-positive <i>S</i> . staphylococci. The y clumping factor) and produced by clinical solution, a rapid aggl fibrinogen w IgG with the The reagents used shift	clumping factor, wl plasma, similarly as ognizes a somatic and occal protein A, and <i>aureus</i> and clum ellow latex particles IgG antibodies that <i>S. aureus</i> isolates. E utination of the parti ith the clumping-fact staphylococcal caps ould be at room temp	hich is specific for a the coagulase. The d a capsular antigen, is serves to distinguish ping factor-negative are coated with fibr react with the type 5 By mixing <i>S. aureus</i> cles occurs by the react for ule.	S. aureus, causes a Staphaurex-Plus slid in addition to clumpin between the clumpin re coagulase-negative inogen (a ligand of the capsule, which is ofte with a drop of the late action of	an de ng ng ve he en ex
For each test add one coated latex) in the se	e drop of test-latex in econd circle of the re	a circle-, and one dr action cardboard.	op of control latex (n	ot
Mix a loopful of the should be suspende precipitate in the test	same bacteria in the o d homogeneously in -latex.	drop of the control- and the control-latex	and test-latex. S. aure. and should produce	us a
Results				
Reaction	Strain 1	Strain 2	Strain 3	
Catalase				
Coagulase				
Clumping factor				

Material and	Material and equipment needed:
Experimental Protocols	Differential and selective plates, sterile cotton tips and loops, glass slides, $H_2O_2$ solution, tube plasma, Staphaurex-Plus kit
	Procedures employed: sampling of bacteria of the nose, fractionation of bacteria on plates, biochemical tests for differentiation
Laboratory Rules &	Laboratory rules:
Precautions	- Do not drink or eat in the room.
	<ul> <li>Persons with skin lesions wear gloves (gloves are provided).</li> </ul>
	- Disinfect the table with 70 % ethanol if you spill organisms.
	<ul> <li>Disinfect your hands before leaving the course.</li> </ul>
Goals & Experiences	Identification of bacteria on selective growth media
gained	biochemical characterization of staphylococci
Timing	$1^{st}$ day, 30 min; $2^{nd}$ day 45 min
Reporting	Results, Discussion, Explanations
Questions to be	How can you distinguish between <i>S aureus</i> and <i>S enidermidis</i> ?
answered	Which strain is more nathogenic: <i>S qureus</i> or <i>S enidermidis</i> ?
	How can you determine experimentally if a strain produces hemolysins?
	What are risk factors for staphylococcal infections?
	How do staphylococci spread from person to person?
	What precautions would you suggest to prevent spreading of staphylococci between patients?
	How many out of all students participating in the practicum are carriers of S. <i>aureus</i> ?