Experiment 10

Bacterial genetic exchange: Conjugation of *ina* (*ice nucleation active* gene from *E. coli* into luminous *V. harveyi*

Authors: Nadja Morf

Gioia Sirena Pirmin Scheuber

Tutor: Munti Yuhana (myuhana@botinst.unizh.ch)

Introduction: For goals and experiences gained look at the instructions in the course

handout.

Practical Work: We get 2 tubes with *E. coli* and 1 tube of *V. harveyi* with the genotypes:

E. coli (pJL 1703): tra, ina, Kan, Amp, lux (donor/D) E. coli (pRK 2013): tra, ina, Kan, Amp, lux (helper/H) V. harveyi: tra, ina, Kan, Amp, lux (recipient/R)

We fill 3 Eppendorf tubes with 1ml D, 3 with 1ml H and 4 with 1ml R each. Then we put them 2 min. in the centrifuge at 7000 rpm and pour off the

supernatant.

4 treatments are controls; R, R+D, R+H and D+H each mixed with 100 μ l *LB*, transfered into 1ml *LB*-agar, incubated overnight at 30°C and then spread on

agar plates (LB, Kan, Amp), next incubation for 48 hours at 30°C.

We give 100 μ I LB into the tube with R mix it and put the liquid into the tube with D and then into the tube with H. At the end we transfer it into a tube with 1 mI LB-medium and incubate it overnight at 30°C. After this we dilute with 0,9 mI LB medium containing Kan and Amp and take 25μ I, 50μ I and 75μ I to spread out on agar plates (LB, Kan, Amp), next incubation for 48 hours at 30°C. In the dark we mark some luminescent colonies and put them on new agar plates (LB, Kan, Amp) that we get enough cells for the ice nucleation assay, for which we take cells and transfer them into PBS. We put the test-

tubes in a cooling bath for 5 min. And look if it freezes.

Results/Discussion: Nothing grew on the control agar plates: R is sensitive to Kan, R+D make no

conjugation because D has *tra-* and D is sensitive to *Amp*, R+H make conjugation but the plasmid isn't working in *V. harveyi*, D+H make

conjugation but they are both sensitive to Amp.

We see luminescence on the agar plates and also the freeze test is positive. So we can say that there is conjugation in some cells in the test-mixture (R+D+H), where the plasmid of H goes into D and helps that both plasmids can get into R. These new cells are resistant to *Kan* and *Amp* and have *ina+*.

Appendix: http://www.microeco.unizh.ch/uni/kurs/bio3_04/start/framepg1-css-b.html

and

http://www.microeco.unizh.ch/uni/kurs/bio3 04/docs/ex-10.html

for the experimental protocols

Reading Chapters 9.1, 9.5, 9.9 in BBOM 9th